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1.1 Iekļautie modeļi

Seller model	Description	Main difference
IRS-HBOT-01	Hard-shell lying type 2.0 ATA chamber - smaller version	Smaller chamber width and lying platform
IRS-HBOT-02	Hard-shell lying type 2.0 ATA chamber - larger version	Larger chamber width and lying platform

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Parameter	IRS-HBOT-01	IRS-HBOT-02 / larger version	Notes
External dimensions	77 x 82 x 225 cm	77 x 90 x 230 cm	Seller model data.
Lying platform dimensions	57 x 220 cm	74 x 225 cm	Seller model data.
Pressure setting range	1.3 / 1.5 / 1.7 / 2.0 ATA	1.3 / 1.5 / 1.7 / 2.0 ATA	Adjustable levels without valve replacement.
Maximum pressure	2.0 ATA absolute / approx. 100 kPa gauge / 14.5 psi gauge	Same	Gauge pressure is above atmospheric pressure.
Chamber material	Composite material	Composite material	Pressure-bearing body.
Hatch / opening size	Approx. 830 x 680 mm where applicable	Approx. 830 x 680 mm where applicable	May vary by model.
Noise level	<= 60 dB(A)	<= 60 dB(A)	From chamber manual.
Pressurisation / depressurisation	Automatic and manual	Automatic and manual	Emergency release available.
Chamber net weight	Approx. 130 kg to be confirmed	Approx. 190 kg to be confirmed	Confirm before installation.
Gross packed weight	Approx. 230 kg to be confirmed	Approx. 290 kg to be confirmed	Confirm before shipment.
Packed chamber size	Approx. 245 x 100 x 105 cm	Approx. 245 x 123 x 117 cm	Packing information.

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Parameter	Value
Oxygen flow rate	Up to 10 L/min, adjustable
Air flow rate	Approx. 110 L/min
Oxygen outlet pressure	>= 350 kPa
Air pressure	>= 200 kPa
Oxygen concentration shown on display	Example display: 93%
Voltage	AC 110/220 V, 50/60 Hz
Power consumption	Approx. 2000 W
Unit dimensions	Approx. 400 x 470 x 1070 mm
Net / gross weight	Approx. 88-89 kg / 100 kg
Packed size	Approx. 60 x 54 x 127 cm
Operating ambient temperature	5 C to 40 C
Operating relative humidity	30% to 75%
Operating atmospheric pressure	860 hPa to 1060 hPa
Display functions	Chamber pressure, oxygen concentration, set pressure, set time, running time, cumulative working time, temperature, humidity, pressurisation and depressurisation rate

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Parameter	Value / note
Function	Cooling water circulation for chamber comfort / air-conditioning circuit
Temperature range	2 C to 10 C; confirm final setting range before commissioning
Input power	Approx. 0.6 kW
Voltage	220/110 V, 50/60 Hz
Unit size	Approx. 365 x 445 x 600 mm
Net / gross weight	Approx. 30 kg / 37 kg
Packed size	Approx. 50 x 48 x 71 cm

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- ATA ir absolūtais spiediens. 2.0 ATA nozīmē aptuveni divkārtšu atmosfēras spiedienu.
- kPa / bar var nozīmēt manometrisko spiedienu virs atmosfēras spiediena.
- 100 kPa gauge ir aptuveni 1.0 bar gauge un atbilst aptuveni 2.0 ATA absolute.
- 14.5 psi gauge ir aptuveni 100 kPa gauge.
- 0.010-0.025 MPa/min ir aptuveni 10-25 kPa/min.

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

External chamber elements

Internal chamber elements

No.	Component	Function
1	Hatch	Entry and exit; must be closed and sealed before pressurisation.
2	Push-type manual pressure relief valve	Manual emergency depressurisation.
3	Pressure gauge	Visual chamber pressure indication.
4	Screw-type manual pressure relief valve	Manual pressure release / adjustment.
5	Handle	Hatch / handling point.
6	Emergency button	Emergency stop / emergency function.
7	Intercom	Communication between user and operator.
8	Internal display	Internal touchscreen / controls.
9	Rail decorative cover	Interior trim component.
10	Chamber air circulator	Internal air circulation.
11	Chamber mattress	User lying platform.
12	Chamber pillow	User comfort accessory.

7.2 Ģeneratora / koncentratora uzbūve

Oxygen generator / concentrator appearance drawing

No.	Component
1	Top cover
2	Screen
3	Shell
4	Caster
5	Switch
6	Flowmeter
7	Overload protection
8	Back cover
9	Air filter cover
10	Aviation socket
11	Pressure detection port
12	Door pneumatic outlet
13	Oxygen outlet
14	Power cord tie
15	Power cord
16	Air outlet
17	Balanced return port
18	Intercom phone

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System flow schematic

8. Pieslēguma punkti


Mark	Connection	Function / corrected description
A	Chamber pressure monitoring interface	Pressure sensing / monitoring connection.
B	Hatch pressurisation / sealing interface	Hatch pressure / pneumatic sealing function.
C	Oxygen input interface	Oxygen supply for approved breathing interface.
D	Chamber air return / back-flow interface	Air return / circulation connection.
E	Air input interface	Main compressed air input to chamber.
1	Oxygen generator connection harness	Electrical / control connection to generator.
2	Cooling water inlet	Cooling unit water inlet connection.
3	Cooling water outlet	Cooling unit water outlet connection.

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


Main touchscreen display example

Screen / symbol	Meaning	Use
Setting key	Language and system settings	Enter settings interface.
Lock screen key	Lock / unlock	Lock or unlock chamber screen.
On / Off key	Start / stop host	Turn host on or off.
Air circulation	Air circulation	Turn internal circulation on/off.
Disinfection key	UV disinfection	Must not be used with a person inside.
Light key	Lighting	Turn chamber light on/off.
Pressure display	Chamber pressure	Real-time pressure value.
O2 purity display	Oxygen concentration	Real-time oxygen concentration display.
Time setting	Session duration	Set session running time.
Pressure setting	Pressure setpoint	Set chamber pressure.
Pressurisation rate	Pressure rise speed	Set pressure increase speed.
Depressurisation rate	Pressure release speed	Set pressure release speed.
Running time	Current session time	Display current session running time.
Cumulative working time	Total working hours	Display cumulative running time.
Humidity / temperature	Environment display	Display chamber humidity and temperature.
Pressure curve	Pressure trend	Display pressure change curve.

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11.5 Maintenance log

Date	Operator	Working hours	Inspection / cleaning performed	Faults / action

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12.1 Chamber troubleshooting

Fault	Possible cause	Action
Chamber does not pressurise	Relief valve open; hatch seal not closed; leakage; generator not started.	Close valve, check hatch, check connections, restart generator.
Louder noise	Uneven placement, vibration or component fault.	Place on level surface; stop use if noise remains.
No oxygen through mask / tube	Loose joint or bent tube.	Insert connector correctly and straighten tube.
Automatic valves abnormal	Control failure or valve damage.	Stop use and contact supplier.
Cannot equalise ears	Pressure rise too fast or ear/sinus issue.	Pause, slow pressurisation or end session.

12.2 Generator alarm codes

Alarm code	Alarm item	Trigger condition	Action
E1	Low chamber pressure	Abnormal low chamber pressure.	Check connections, hatch seal, valves and tubing.
E2	High chamber pressure	Abnormal high chamber pressure.	Stop session, depressurise safely and contact supplier if repeated.
E3	Oxygen system failure	Oxygen concentration is low.	Check oxygen flow, filters and vents. Contact supplier if repeated.
E4	Emergency stop fault	Emergency stop is pressed.	Reset only after confirming safe conditions.
E5	Communication failure	Signal interruption.	Stop use if control communication is not restored.

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15.1 Chamber accessories

Accessory	Typical quantity	Notes
Mattress	1 set	For chamber interior.
Pillow	1 set	For user comfort.
Intercom	1 pc	Communication system.
Oxygen mask	5 pcs	Quantity may vary by order.
Headset tube / nasal oxygen tube	1 pc	Quantity may vary by order.
User manual	1 pc	This document or current approved version.

15.2 Generator accessories

Accessory	Typical quantity	Notes
Main generator / concentrator	1 pc per system	Quantity may vary.
Primary air intake filter	1 pc	Clean weekly.
Secondary air intake filter	2 pcs	Replace at least every 3,000 hours.
User manual / certification	1 set	Internal technical file / supplied documents.

15.3 Packed component weights

Component	Approx. net weight	Approx. gross weight	Approx. packing size
Oxygen generator / concentrator	88-89 kg	100 kg	60 x 54 x 127 cm
IRS-HBOT-01 chamber	130 kg	230 kg	245 x 100 x 105 cm
IRS-HBOT-02 chamber	190 kg	290 kg	245 x 123 x 117 cm
Cooling water circulation unit	30 kg	37 kg	50 x 48 x 71 cm

16. Operatora kontrolsaraksts

16.1 Pirms seansa

Check	OK / Not OK	Notes
Safety restrictions explained to user.		
No obvious no-use condition or unresolved health concern.		
Prohibited items removed.		
Cotton or antistatic clothing confirmed.		
Chamber, hatch, window, valves and seals inspected.		
Generator, cooling unit, tubing and connections inspected.		
Intercom / communication verified.		
Pressure, time and rate settings checked.		
Operator remains available for full session.		

16.2 Pēc seansa

Check	OK / Not OK	Notes
User exited safely; no unusual discomfort reported.		
Pressure returned to zero before hatch opening.		
Power switches and valves set correctly after use.		
Chamber cleaned and disinfected as required.		
Breathing accessories cleaned/replaced according to hygiene rules.		
Faults, alarms or abnormal behaviour recorded.		

Dokumenta kontrole

Item	Value
Document	IRS-HBOT cietā korpusa 2.0 ATA piņas sistēmas lietotāja rokasgrāmata
Language	Latviešu
Use type	Non-medical wellness use only
Document owner	UAB OMETA
Revision	LV v0.4 tehniskai pārbaudei
Note	Final release requires technical confirmation of exact model specifications, PED/pressure equipment classification and final packing list.